



Clean Pathogen DNA & RNA Kit

Instructions for Use

V.9 - MAY 2024

For Research Use Only

REF CPT-DR0096, CPT-DR0384

 CleanNA, Coenecoop 75, 2741 PH, Waddinxveen, The Netherlands

Information in this document is subject to change without notice.

Disclaimer

For Research Use Only. Not for use in diagnostic procedures.

CleanNA disclaims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. To the fullest extent allowed by law, in no event shall CleanNA be liable, whether contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including, but not limited to the use thereof, whether or not foreseeable and whether or not CleanNA is advised of the possibility of such damages.

Trademarks

The trademarks mentioned herein are the property of their respective owners.

Contents

Intended Use	4
Introduction and Principle	4
Schematic Overview	5
Materials Provided.....	6
Reagent Shipping, Storage and Handling	6
Warnings	7
Precautions.....	8
Quality Control	9
Materials and Equipment to be Supplied by User	10
Working RNase Free	10
Preparation of Reagents	11
Tissue Protocol	12
Serum & Stool Protocol.....	15
Troubleshooting Guide	18
Symbols	20
Ordering Information	21
Document Revision History	22
Notes	23

Intended Use

Clean Pathogen DNA/RNA Kit is intended for use by professional users trained in molecular biology techniques. It is designed to use manually or on a liquid handling workstation for molecular biology applications.

Introduction and Principle

The Clean Pathogen DNA/RNA Kit is designed for high throughput and reliable isolation of both DNA and RNA from tissue, serum and fecal samples.

The Clean Pathogen Kit will isolate DNA from:

- High quality host genomic DNA
- Gram positive and negative bacterial DNA
- Fungal spore DNA
- Viral DNA

The Clean Pathogen Kit will isolate RNA from:

- Host RNA
- Viral RNA

The Clean Pathogen DNA/RNA Kit combines our propriety buffer system with the convenience of our magnetic CleanNA Particles XP. Our buffer system minimizes the binding of PCR inhibiting compounds, present within the samples, onto our magnetic particles.

The protocol is fully scalable and due to the use of our magnetic bead purification technology, can besides manual usage, easily be automated once the samples have been lysed on our CleanXtract 96 or liquid handling workstations (e.g. Dynamic Devices LYNX™, Hamilton STAR™).

Purified DNA is suitable for PCR, qPCR, restriction digestion, Next Generation Sequencing and hybridization applications. Purified RNA is suitable for reverse transcription for cDNA synthesis, RT-PCR, RT-qPCR and RNA-Seq.

Schematic Overview

With our specially formulated lysis buffer in combination with the glass bead mix in the included disruptor plate, the samples are lysed. The DNA and RNA is bound to our magnetic particles while DNases and RNases are deactivated. DNA and RNA is isolated from the lysates in one step by binding to the CleanNA Particles' surfaces. The CleanNA magnetic particles are separated from the lysates by using a magnetic separation device. Following a few rapid wash steps to remove trace contaminants, the purified DNA/RNA is eluted from the CleanNA particles for downstream applications using an Elution Buffer.




Figure 1: Schematic overview of the Clean Pathogen DNA & RNA Kit extraction procedure.

Materials Provided


Kit Contents:

Component	CPT-DR0096 (1x96 preps)	CPT-DR0384 (384 preps)
Clean Disruptor Plate 96	1	4
Bag of caps	1	4
TLA Lysis	60 mL	240 mL
NLB Buffer	8 mL	30 mL
Proteinase K (20 mg/mL)	2.2 mL	9 mL
YB Buffer	40 mL	160 mL
BA Binding	40 mL	160 mL
CleanNA Particles XP	2.2 mL	9 mL
AR Wash	88 mL	3 x 88 mL
EX Wash	30 mL	4 x 30 mL
Elution Buffer	15 mL	60 mL

Reagent Shipping, Storage and Handling

Clean Pathogen DNA & RNA Kit is shipped at room temperature (15-25 °C). Do not freeze the components of the Clean Pathogen DNA & RNA Kit. After the components have been frozen, the kit is no longer suitable for use. Do not use the Clean Pathogen DNA & RNA Kit after the expiration date stated on the outer box label.

Component	Storage Temperature
Clean Disruptor Plate 96	15-25 °C
Caps	15-25 °C
TLA Lysis	15-25 °C
NLB Buffer	15-25 °C
Proteinase K (20 mg/mL)	15-25 °C (for storage > 12 months, store at 2-8 °C)
YB Buffer	15-25 °C
BA Binding	15-25 °C
CleanNA Particles XP	2-8 °C
AR Wash	15-25 °C
EX Wash	15-25 °C
Elution Buffer	15-25 °C

 **Note:** Check all buffers for precipitates prior to usage. Any precipitates can be re-dissolved by warming the buffer(s) to 37°C and shaking gently.

Warnings

Read the instructions carefully before using the kit.

Do not mix several kit LOT numbers.

The LOT number on the CleanNA Particles XP box packaging is different from the LOT number on the CleanNA Particles XP bottle. The LOT number on the box matches the LOT number of the whole kit and the one on the bottle is specifically for the particles. Since the CleanNA Particles XP are stored at a different temperature than the rest of the kit, please make sure that the LOT number on the box packaging of the particles matches the LOT number of the kit before use.

Precautions

When working with chemicals, always follow your facility's procedures and universal precautions by using disposable gloves, safety glasses, a labcoat etc.

For all safety information, please consult the safety data sheet (SDS).

NLB Buffer



Causes serious eye damage. Causes skin irritation.
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor/physician/first aider.
IF ON SKIN: Wash with plenty of water and soap.
If skin irritation occurs: Get medical advice/attention.
Take off contaminated clothing and wash it before reuse.

Proteinase K (20 mg/mL)



May cause an allergic skin reaction. Causes skin irritation. Causes serious eye irritation. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause respiratory irritation.
Wear protective gloves/protective clothing/eye protection/face protection.
Avoid breathing dust/fume/gas/mist/vapours/spray.
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.

YB Buffer



Harmful if swallowed. Causes severe skin burns and eye damage. Harmful to aquatic life with long lasting effects.
IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor/physician/first aider. Wash contaminated clothing before reuse.
IF SWALLOWED: Call a POISON CENTER/doctor/physician/first aider if you feel unwell.
IF INHALED: Remove person to fresh air and keep comfortable for breathing.

BA Binding



Toxic to aquatic life with long lasting effects. Causes serious eye damage. Harmful if swallowed. Causes skin irritation.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor/physician/first aider. Collect spillage.

IF SWALLOWED: Call a POISON CENTER/doctor/physician/first aider if you feel unwell.

IF ON SKIN: Wash with plenty of water and soap. Rinse mouth. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.



AR Wash



Harmful if swallowed. Causes skin irritation. Causes serious eye irritation.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

IF SWALLOWED: Call a POISON CENTER/doctor/physician/first aider if you feel unwell.

IF ON SKIN: Wash with plenty of water and soap. Rinse mouth. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

Note: For safe disposal, please consult your local waste regulations.

Quality Control

CleanNA produces each lot of Clean Pathogen DNA & RNA Kit according to predetermined and validated protocols in the Quality Management System (QMS). CleanNA's QMS is EN-ISO 13485 certified.

Materials and Equipment to be Supplied by User

Materials and reagents to be supplied by user for Tissue and Serum & Stool protocol:

- Centrifuge capable of at least 3,500 x g with adaptor for 96-well plates
- Magnetic separation device for 96-well plates (CleanNA, Part# CMAG-96-RN50)
- Incubator capable of 70°C
- 96-well plates with a capacity of at least 1.7 mL (Recommend Nunc 278752) and compatible with the Magnetic Separation Device
- 96-well microplates for DNA storage
- Vortexer
- Absolute ethanol
- Molecular biology grade water
- Optional: Mixer mill such as a SPEX CertiPrep Geno/Grinder® 2010 or Qiagen TissueLyser

Working RNase Free

For RNA applications it is important to work RNase free. RNases are present everywhere and general precautions should be taken to avoid the introduction of RNases and other contaminating nucleases while working with RNA. The most common sources of RNases are hands, dust particles and contaminated laboratory solutions, equipment and glassware.

To minimize the risk of RNase contamination we recommend the following precautions:

- Always use gloves when handling RNA samples. Change your gloves frequently, to avoid contaminations;
- Ensure to use RNase free filter tips for pipetting;
- Use materials such as disposable consumables, which are guaranteed RNase free;
- Use reagents which are guaranteed RNase free. Creating aliquots from buffers lowers the risk of RNase contamination in buffers, reagents, etc.;
- Avoid using reagents, consumables and equipment dedicated for common use or general lab processes;
- If possible work in a separate room, fume hood or lab space;
- Clean all working surfaces with commercial RNase inhibiting surfactant or 70% ethanol before starting your work / experiment.

Preparation of Reagents

AR Wash

Prepare AR Wash with absolute ethanol as follows and store at room temperature.

Kit	Absolute ethanol to be added
CPT-DR0096	112 mL
CPT-DR0384	112 mL per bottle

EX Wash

Dilute EX Wash with absolute ethanol as follows and store at room temperature.

Kit	Absolute ethanol to be added
CPT-DR0096	70 mL
CPT-DR0384	70 mL per bottle

Shake or vortex the CleanNA Particles XP to fully resuspend the particles prior to usage. The particles must be fully suspended during use to ensure proper binding.

Tissue Protocol


This protocol describes the procedure for the isolation of both host and pathogen DNA & RNA from tissue samples.

Before Starting:


- Prepare AR Wash and EX Wash according to the “Preparation of Reagents” section on Page 11.
- Set an incubator to 70°C.
- Heat Elution Buffer to 70°C.

Protocol:


1. Briefly spin the Clean Disruptor Plate 96 to remove any glass beads from the walls of the wells. Uncap the Clean Disruptor Plate 96 and save the caps for use in Step 3.
2. Add 25-30 mg tissue to each well.
3. Add 525 μ L TLA Lysis to each sample. Seal the Clean Disruptor Plate 96 with the caps removed in Step 1.
4. Vortex at maximum speed for 3-5 minutes to lyse and homogenize samples. For best results, a Mixer Mill, such as Spex CertiPrep Geno/Grinder® 2010 or Qiagen Tissuelyser, should be used.

 **Note:** Depending on the sample amount and type, the amount of TLA Lysis may need to be adjusted so that 300 μ L can be recovered after Step 11.


5. Centrifuge at 1,000-2,000 x g for 60 seconds at room temperature.
6. Remove and discard the caps from the Clean Disruptor Plate 96.
7. Add 53 μ L NLB Buffer and 20 μ L Proteinase K (20 mg/mL) to each sample.
8. Seal the Clean Disruptor Plate 96 with new Caps for Racked Microtubes (provided).
9. Vortex for 60 seconds to mix thoroughly.
10. Incubate at 70°C for 15 minutes. Mix once during incubation.
11. Centrifuge at 3,500 x g for 10 minutes.
12. Transfer 300 μ L cleared supernatant to a 96-well deep-well plate (1.7 mL compatible with the Magnetic Separation Device used).

 **Note:** Do not transfer any debris as it can reduce yield and purity.


13. Add 300 μ L YB Buffer, 300 μ L BA Binding, and 20 μ L CleanNA Particles XP to each sample. Vortex to mix thoroughly or pipet up and down 20 times.

 **Note:** CleanNA Particles XP and YB Buffer can be prepared as a master mix prior to use. Prepare only what is needed. Tip mixing is recommended for automated protocols for best yield.

14. Incubate at room temperature for 10 minutes.
15. Place the 96-well deep-well plate on the Magnetic Separation Device to magnetize the CleanNA Particles XP. Incubate at room temperature until the CleanNA Particles XP are completely cleared from solution.
16. Aspirate and discard the cleared supernatant. Do not disturb the CleanNA Particles XP.
17. Remove the plate containing the CleanNA Particles XP from the Magnetic Separation Device.
18. Add 600 μ L AR Wash to each sample. Resuspend the CleanNA Particles XP by vortexing or pipetting up and down 20 times.

 **Note:** AR Wash must be diluted with absolute ethanol prior to use. Please see Page 11 for instructions.


19. Incubate at room temperature for 2 minutes.
20. Place the 96-well deep-well plate on the Magnetic Separation Device to magnetize the CleanNA Particles XP. Incubate at room temperature until the CleanNA Particles XP are completely cleared from solution.
21. Aspirate and discard the cleared supernatant. Do not disturb the CleanNA Particles XP.
22. Remove the plate containing the CleanNA Particles XP from the Magnetic Separation Device.
23. Repeat Steps 18-22 once for a second AR Wash step.
24. Add 600 μ L EX Wash to each sample. Resuspend the CleanNA Particles XP by vortexing or pipetting up and down 20 times.

 **Note:** AR Wash must be diluted with absolute ethanol prior to use. Please see Page 11 for instructions.

25. Incubate at room temperature for 2 minutes.
26. Place the 96-well plate on the Magnetic Separation Device to magnetize the CleanNA Particles XP. Incubate at room temperature until the CleanNA Particles XP are completely cleared from solution.

27. Aspirate and discard the cleared supernatant. Do not disturb the CleanNA Particles XP.
28. Leave the plate on the Magnetic Separation Device.
29. Add 500 µL molecular biology grade water (not provided) to each sample. Immediately aspirate and remove the molecular biology grade water. Do not let the samples stay in contact with the molecular biology grade water for more than 60 seconds.

Note: This step (water “wash”) is crucial to remove traces of ethanol from the sample. When adding water, the CleanNA Particles XP will come out of the ring shape shortly. This is a normal phenomenon, the CleanNA particles XP will return to the ring shape after a few seconds.

 **Note:** If using an automated platform, use the maximum volume the tips will allow up to 600 µL.

30. Add 50-100 µL Elution Buffer heated to 70°C to each sample. Resuspend CleanNA Particles XP by vortexing or pipetting up and down 20 times.
31. Incubate at room temperature for 5 minutes.
32. Place the 96-well plate on the Magnetic Separation Device to magnetize the CleanNA Particles XP. Incubate at room temperature until the CleanNA Particles XP are completely cleared from solution.
33. Transfer the cleared supernatant containing purified DNA/RNA to a clean 96-well microplate. Store the DNA at -20°C. For downstream RNA applications, store the eluate at -80°C.

Serum & Stool Protocol

This protocol describes the procedure for the isolation of both host and pathogen DNA & RNA from serum and stool samples.

Before Starting:

- Prepare AR Wash and EX Wash according to the “Preparation of Reagents” section on Page 11.
- Set an incubator to 70°C.
- Heat Elution Buffer to 70°C.

Protocol:

1. Briefly spin the Clean Disruptor Plate 96 to remove any glass beads from the walls of the wells. Uncap the Clean Disruptor Plate 96 and save the caps for use in Step 3.
2. Add 250 μL serum or stool samples to each well. If stool sample is solid, resuspend to 10% wgt/volume in PBS before starting.
3. Add 275 μL TLA Lysis to each sample. Seal the Clean Disruptor Plate 96 with the caps removed in Step 1.
4. Vortex at maximum speed for 3-5 minutes to lyse and homogenize samples. For best results, a Mixer Mill, such as Spex CertiPrep Geno/Grinder® 2010 or Qiagen Tissuelyser, should be used.



Note: Depending on the sample amount and type, the amount of TLA Lysis may need to be adjusted so that 300 μL can be recovered after Step 11.

5. Centrifuge at 1,000-2,000 x g for 60 seconds at room temperature.
6. Remove and discard the caps from the Clean Disruptor Plate 96.
7. Add 50 μL NLB Buffer and 20 μL Proteinase K (20 mg/mL) to each sample.
8. Seal the Clean Disruptor Plate 96 with new Caps for Racked Microtubes (provided).
9. Vortex for 60 seconds to mix thoroughly.
10. Incubate at 70°C for 15 minutes. Mix once during incubation.
11. Centrifuge at 3,500 x g for 10 minutes.

- Transfer 300 μ L cleared supernatant to a 96-well deep-well plate (1.7 mL) compatible with the Magnetic Separation Device used.



Note: Do not transfer any debris as it can reduce yield and purity.

- Add 300 μ L YB Buffer, 300 μ L BA Binding, and 20 μ L CleanNA Particles XP to each sample. Vortex to mix thoroughly or pipet up and down 20 times.



Note: CleanNA Particles XP and YB Buffer can be prepared as a master mix prior to use. Prepare only what is needed. Tip mixing is recommended for automated protocols for best yield.

- Incubate at room temperature for 10 minutes.
- Place the 96-well deep-well plate on the Magnetic Separation Device to magnetize the CleanNA Particles XP. Incubate at room temperature until the CleanNA Particles XP are completely cleared from solution.
- Aspirate and discard the cleared supernatant. Do not disturb the CleanNA Particles XP.
- Remove the plate containing the CleanNA Particles XP from the Magnetic Separation Device.
- Add 600 μ L AR Wash to each sample. Resuspend the CleanNA Particles XP by vortexing or pipetting up and down 20 times.



Note: AR Wash must be diluted with absolute ethanol prior to use. Please see Page 11 for instructions.


- Incubate at room temperature for 2 minutes.
- Place the 96-well deep-well plate on the Magnetic Separation Device to magnetize the CleanNA Particles XP. Incubate at room temperature until the CleanNA Particles XP are completely cleared from solution.
- Aspirate and discard the cleared supernatant. Do not disturb the CleanNA Particles XP.
- Remove the plate containing the CleanNA Particles XP from the Magnetic Separation Device.
- Repeat Steps 18-22 once for a second AR Wash step.
- Add 600 μ L EX Wash to each sample. Resuspend the CleanNA Particles XP by vortexing or pipetting up and down 20 times.



Note: EX Wash must be diluted with absolute ethanol prior to use. Please see Page 11 for instructions.

- Incubate at room temperature for 2 minutes.

26. Place the 96-well plate on the Magnetic Separation Device to magnetize the CleanNA Particles XP. Incubate at room temperature until the CleanNA Particles XP are completely cleared from solution.
27. Aspirate and discard the cleared supernatant. Do not disturb the CleanNA Particles XP.
28. Leave the plate on the Magnetic Separation Device.
29. Add 500 μ L molecular biology grade water (not provided) to each sample. Immediately aspirate and remove the molecular biology grade water. Do not let the samples stay in contact with the molecular biology grade water for more than 60 seconds.

 **Note:** This step (water “wash”) is crucial to remove traces of ethanol from the sample. When adding water, the CleanNA Particles XP will come out of the ring shape shortly. This is a normal phenomenon, the CleanNA particles XP will return to the ring shape after a few seconds.

Note: If using an automated platform, use the maximum volume the tips will allow up to 600 μ L.

30. Add 50-100 μ L Elution Buffer heated to 70°C to each sample. Resuspend the CleanNA Particles XP by vortexing or pipetting up and down 20 times.
31. Incubate at room temperature for 5 minutes.
32. Place the 96-well plate on the Magnetic Separation Device to magnetize the CleanNA Particles XP. Incubate at room temperature until the CleanNA Particles XP are completely cleared from solution.
33. Transfer the cleared supernatant containing purified DNA/RNA to a clean 96-well microplate. Store the DNA at -20°C. For downstream RNA applications, store the eluate at -80°C.

Troubleshooting Guide







Please use this guide to troubleshoot any problems that may arise. For further assistance, please contact your local distributor.

Possible problems and Suggestions

Problem	Cause	Suggestion
Low DNA yield or no DNA yield	Poor homogenization of sample.	Repeat the DNA isolation with a new sample, be sure to mix the sample with TLA Lysis thoroughly. Use a commercial homogenizer if possible.
	Incomplete resuspension of CleanNA Particles XP.	Resuspend the CleanNA Particles XP by vortexing vigorously before use.
	DNA washed off.	Make sure AR Wash and EX Wash are mixed with ethanol.
		Ensure the water "wash" does not exceed 60 seconds and the CleanNA Particles XP are not resuspended.
CleanNA Particles XP lost in process.	After water is added during wash step CleanNA Particles XP will go into solution. Let magnetic particles re-magnetize prior to aspirating liquid.	
cDNA synthesis not working	Degraded RNA.	Store the eluate at -80°C after isolation and minimize the number of freeze thaw cycles.
		Ensure to work RNase free and use RNase free materials during the isolation process and in downstream applications.
CleanNA Particles XP do not completely clear from solution	Magnetizing time too short.	Increase the particle collection time on the magnetic separation device.
A260/A230 ratio is low	Salt contamination.	Repeat the DNA isolation with a new sample. Extend the incubation time with AR Wash. Wash the CleanNA Particles XP with ethanol.
A260/A280 ratio is high	RNA contamination.	The protocol does not remove RNA. If desired, add 5 µL RNase A (25 mg/mL) (not provided) after lysate is cleared and before binding buffers are added. Incubate at room temperature for 5 minutes.

Problem	Cause	Suggestion
Problems in downstream applications	BSA not added to PCR mixture.	Add BSA to a final concentration of 0.1 µg/mL to the PCR mixture.
	Too much DNA inhibits PCR reactions.	Dilute the DNA elute used in the downstream application if possible.
	Non-specific bands in downstream PCR.	Use hot-start Taq polymerase mixture.
	Inhibitory substance in the eluted DNA.	Check the A260/A230 ratio. Dilute the elute to 1:50 if necessary.
Abnormal BioAnalyzer data	BioAnalyzer shows multiple sharp peaks during the analysis.	Store the eluate at -80°C after isolation and minimize the number of freeze thaw cycles.
		Ensure to work RNase free and use RNase free materials during the isolation process and in downstream applications.
	BioAnalyzer shows base line climbing towards the end.	Check the BioAnalyzer chip for air bubbles. Load samples onto a new freshly prepared chip.
	BioAnalyzer shows high blob at the beginning of the trace.	Ensure the purified sample does not contain traces of CleanNA Particles XP.

Symbols

	Reference number
	Manufacturer
	Caution
	Temperature limit
	Expiration date
	Lot number

Ordering Information

Contact your local distributor to order.

Product	Part Number
Clean Pathogen DNA & RNA Kit 96 preps	CPT-DR0096
Clean Pathogen DNA & RNA Kit 384 preps	CPT-DR0384

Product	Part Number
Clean Magnet Plate 96-well RN50	CMAG-96-RN50

Document Revision History

Manual Version	Date of revision	Revised Chapter	Explanation of revision
9	07/MAY/2024	Total Revision	Updated buffer names and volumes.
8	29/FEB/2024	Materials Provided and Urine and Whole Blood protocol	Removed Urine and Whole Blood protocol from manual.
7	27/OCT/2023	Total revision.	New lay-out. Added chapters schematic overview, warnings, precautions, quality control, symbols.
6.00	October 2021	Total revision.	Language and layout modifications.
		Kit contents and material.	Added CPT+ Reagent to table.
5.00	August 2020	Total revision.	New lay-out.
			Important general information added at page 1 (before contents).

Notes

Contact

Coenecoop 75 | 2741 PH Waddinxveen | The Netherlands

T: +31 (0) 182 22 33 50 | F: +31 (0) 182 22 33 98 | info@cleanna.com

www.cleanna.com

